



2023 PATHOLOGY PRODUCT CATALOGUE

Guangzhou Wexis Biotech Limited

No.4 of Yin Sha Road 2,Xin Tang,Zeng Cheng
District,Guangzhou,China
Tel: +86 20 87793718
Fax: +86 20 87793728
Web: www.wexisgp.com
E-mail: marketing@wexisgp.com

Wexis Group Limited

Unit1611,16/FI,Yuen Long Trading Centre,33
Wang Yip Street,Yuen Long,NT,Hong Kong
Tel: +852 25672882
Fax: +852 25676333



*Images are for illustration purposes only
CY-GHML-2301E



2023 PATHOLOGY PRODUCT CATALOGUE

Guangzhou Wexis Biotech Limited

COMPANY PROFILE

Guangzhou Wexis Biotech Limited founded in 2005, integrates R & D, production, sales and service into one, established professional pathological products and overall solutions, and layout nationwide, radiation overseas. After years of accumulation and precipitation, the company has become a national high-tech enterprise in the production and research of conventional pathological products at home and abroad.

We have passed the ISO13485 Medical Device Quality Management System Certification, Intellectual Property Management System Certification, Two Integration Management System (Grade A) Certification, Guangzhou Cleaner Production Enterprise Certification, and in 2017 became the "Pathology Special Training Base" designated by the Pathology Equipment Branch of China Medical Equipment Association.

CONTENTS

01	CRYOTOME
03	ROTARY MICROTOME
05	DRYING CABINET
07	TISSUE SECTION ATOMIZER
09	COLD PLATE
11	WAX TRIMMER
13	FLOATATION / DRYING WORKSTATION
15	SLIDE DRYING HOTPLATE
17	WAX DISPENSER
19	PATHOLOGY SPECIMEN FILING SYSTEM
21	INTENSIVE PATHOLOGICAL SPECIMEN FILING CABINET
23	GROSSING STATION
25	CRYO-CABINET

01

27	EMBEDDING CASSETTES
31	STAINLESS STEEL EMBEDDING MOLDS
31	STAINLESS LIDS
33	DRYING PLATE
35	PATHOLOGY GRADE PARAFFIN WAX
37	DISSECTION BLADE HANDLE
39	MICROSCOPE SLIDES
41	COVER GLASS
42	FROZEN TISSUE EMBEDDING CENTER

02

43	TISSUE FIXATIVE
45	HEMATOXYLIN-EOSIN SOLUTION
47	SPECIAL STAIN SOLUTION
49	HYALINE LIQUID
49	DEHYDRATED LIQUID
50	OTHER REAGNT

03

CRYOTOME C5500

- Dual compressor refrigeration work, accurate temperature control.
- Handwheel running smoothly, can be locked when it be needed.
- Electric heated window can effectively prevent fogging.
- Spot UV light source provides effective disinfection function.
- HD color LCD touch screen.
- Equipped with debris tray and shelving.
- Quick freezing bar can hold 24 specimens at a time.



SPECIFICATION PARAMETERS

Item No.	313070005
Dimensions	760mm*1200mm*835mm
Weight	140kg
Voltage	AC110V~220V/50Hz
Power Rating	1200W
Chamber Temperature	0°C~-60°C
Frosting	Manual/Automated
Cryobar Minimum Temperature	-55°C
Cryobar Freezing Points	≥24pcs
Cryobar Semiconductor Freezing Points	6pcs
Specimen Horizontal Stroke	22mm
Specimen Vertical Stroke	65mm
Max Specimen Size	65mm*80mm
Max Specimen Clamp Adjustable Angle	12°
Specimen Fast Forward And Fast Backward Speed	0.45mm/s
Section Thickness Set Range	0.5μm~100μm
Trim Thickness Set Range	0μm~600μm

ROTARY MICROTOME M3050

- Movable control panel, HD color touch screen.
- Handle is made of antibacterial material.
- Hand wheel lock and 360° any position lock dual locking system.
- Automatically select the count of section from section and trim.
- Highly accurate sample feed system, Precise dynamic balance handwheel design.
- Special process, anti-static, non-stick wax



SPECIFICATION PARAMETERS

Item No.	313070006
Dimensions	580mm*443mm*320mm
Specimen Vertical Stroke	60mm
Specimen Horizontal Stroke	28±1mm
Specimen Positioning	X/Y:0°~10°,Z:360°
Specimen Clamp Adjustable Angle	0°~15°
Section Thickness Set Range	0~100μm, Graduation value 0.1μm
Trim Thickness Set Range	0~999μm, Graduation value 1μm
Specimen Retraction	0~100μm

DRYING CABINET

- Ventilation design, fast drying.
- Orderly display, work more organized.
- Environmentally friendly and healthy.
- Mobile cart, reasonable use a space.

SPECIFICATION PARAMETERS



Item No.	313130001
Product Name	Tabletop Drying Cabinet
Model	PC7010
Configuration	NCCO1804
Capacity	200 Slides
Dimensions	235mm*340mm*290mm



Item No.	313130002
Product Name	Tabletop Drying Cabinet
Model	PC7020
Configuration	NCCO1804
Capacity	400 Slides
Dimensions	235mm*340mm*448mm

SPECIFICATION PARAMETERS

Item No.	313130003	313130004
Product Name	Base type Drying Cabinet	Standard type Drying Cabinet
Model	EC700-1	EC700-2
Configuration	2 Cabinets + Chassis	Base type +80 Drying plates
Capacity	1600 Slides	1600 Slides

Item No.	313130005	313130006
Product Name	Environment Friendly	Environment Friendly Heating
Model	EC700-3	EC700-4
Configuration	Standard type+NCCO1804	Environment Friendly+ Heating
Capacity	1600 Slides	1600 Slides



TISSUE SECTION ATOMIZER NE410

- Cooperate with rotary microtome used.
- Say goodbye to broken, rolled, static.
- Ultrasonic atomization technique.
- Fog volume and air volume can be adjusted.
- Fog tube randomly adjust the fog le of the fuselage.
- Water shortage protection, led water level status indicated.
- Small, light, with no extra space.



SPECIFICATION PARAMETERS

Item No.	313070003
Dimensions	265.5mm*138.5mm*162.6mm
Weight	1.5kg
Atomization Capacity	350ml/h
Voltage	DC24V
Power Rating	0~28W

COLD PLATE CP380

- WEXIS' s CP380 is an ideal assistant in pathology sectioning process.
- Effectively maintaining paraffin block temperature for smooth sectioning.
- Eco-compressor refrigeration system, long service life & efficient freezing.
- Unique groove design to avoid water leaking from cryo-workbench.



SPECIFICATION PARAMETERS

Item No.	313070001
Dimensions	535mm*330mm*358mm
Weight	20kg
Cooling plate area	270mm*250mm
Line-in(voltage)	220V.50Hz;110V.60Hz
Rating	170W
Temperature Range	-2°C~-15°C
Frost Time	≥5 mins
Cooling Plate Capacity	48pcs

WAX TRIMMER TR302

- Low centre of gravity, with ergonomically designed hand rest, operator is able to maintain firm positioning of the unit during operation.
- Simple and pleasing design, the unit is water proof, well insulated and heat resistant.
- Unique heating system, fast heat up with accurate and constant temperature.
- Tilted aluminum alloy heating surface with grooves, offers smooth drainage of excess waste wax.
- Built-in waste collection drawer, for ease and convenient waste collection and disposal.



SPECIFICATION PARAMETERS

Item No.	313070002
Dimensions	285mm*160mm*120mm
Weight	1.95kg
Line-in(voltage)	AC230V~50Hz, 0.3A
Rating	50W
Surface Temperature	78°C (±2°C)
Preheat Time	5~8 mins
The Capacity of Wax Collection Box	300ml
Working Temp	5°C~40°C

FLOATATION/DRYING WORKSTATION FL220

- Combining slide floatation and drying in one body.
- Front panel is tilted 40 degree with IC touch-sensitive panel, provides users with ease and comfort operation.
- Special gill designed PE outer casing, avoid water splashes and ideally insulated, offering safe and relaxed operation.
- PTFE water bath and drying panel are both treated with dark color, ideal contrast offers easy observation of para sections.
- Separate control for oat bath and dryer, LED digital display.
- Detachable PTFE water bath - water replacements.



SPECIFICATION PARAMETERS

Item No.	313080001
Dimensions	415mm*310mm*120mm
Floatation Bath	220mm*145mm
Dryer Panel Area	248mm*107mm
Weight	5.45kg
Line-in(voltage)	AC220V~50Hz 0.3A
Rating	260W(Bath-200W;Dryer-60W)
Dryer Temp Set	Room Temp -99°C
Bath Temp Set	Room Temp -99°C
Working Temp	0°C~40°C
Time of Reach Set Temp	<30 mins

SLIDE DRYING HOTPLATE DR230

- Large hotplate surface.
- Slim and light Design.
- Built-in Micro-processor Controlled Heating System.
- “Stacking Tiles” Style 4 beveled sections hotplate.
- Drying hotplate with large capacity dries a maximum of 44 slides at one time.
- LED digital display, temperature settings from 0°C to 99°C with auto memory of last settings.
- Accurate temperature controls with automatic adjustments according to set temperature.



SPECIFICATION PARAMETERS

Item No.	313080002
Dimensions	415mm*310mm*120mm
Dryer panel area	295mm*270mm
Weight	6.2kg
Line-in(voltage)	AC220V~50Hz 0.3A
Rating	300W
Dryer Temp Set	Room Temp -99°C
Temp Display Range	0°C~99.9°C
Working Temp	0°C~40°C
Time of Reach Set Temp	≤10 mins

WAX DISPENSER WD360

- 6L large-capacity aluminum alloy wax cylinder, melting 4.5Kg of solid wax at a time.
- Microprocessor control heating temperature, continuous heating with accurate constant temperature.
- Heated faucet to ensure smooth flow of wax, without blocking.
- Built-in filter at faucet outlet, minimize impurities outflow.
- Extra heat insulation cover, offering safe operation.



SPECIFICATION PARAMETERS

Item No.	313060001
Dimensions	248mm*289mm*460mm
Weight	9kg
Line-in(voltage)	AC220~50Hz
Wax Cylinder Capacity	6L
Melting Wax Weight	4.5kg
Temperature Setting	0~99°C
Working Temp	0°C~40°C
Wax Outflow Time	<10 mins

PATHOLOGY SPECIMEN FILING SYSTEM



SLIDE STORAGE SYSTEM

Combined by 4 small cabinets;9 drawers per small cabinet,total 36 drawers;each drawer provides 4 rows for glass slide storage.



DRY SLIDE STORAGE SYSTEM

Combined by 4 small cabinets;5 drawers per small cabinet,total 20 drawers;each drawer provides 5 rows for slide drying post medium mounting.



WAX BLOCK FILING SYSTEM

Combined by 4 small cabinets;5 drawers per small cabinet,total 20 drawers; each drawer provides 14 rows.



SPECIFICATION PARAMETERS

Item No.	313130007
Product Name	Slide Storage System
Model	DC700
Dimensions	515mm*480mm*1440mm
Capacity	55000-65000pcs
Loading Weight	1000kg

Item No.	313130008
Product Name	Dry Slide Storage System
Model	HC700
Dimensions	515mm*480mm*1440mm
Capacity	14400-16000pcs
Loading Weight	1000kg

Item No.	313130009
Product Name	Wax Block Filing System
Model	BC700
Dimensions	515mm*480mm*1440mm
Capacity	15000-16500pcs
Loading Weight	800kg

INTENSIVE PATHOLOGICAL SPECIMEN FILING CABINET

- Drawer adopts three steel ball silent slide rails.
- Cabinet surface anti-rust treatment, exquisite appearance.
- Specifies the slot for metal mold opening.
- Integrated stamping forming.
- Centralized management with high usage.



SPECIFICATION PARAMETERS

Product Name	Wax Block Filing System
Dense Frame Dimensions	900mm*500mm*2005mm
Built-in Dimensions	403mm*478mm*1600mm
Handle	Plastic hidden handle
Wax Block Storage Capacity	28800 blocks/group



Product Name	Dry Slide Storage System
Dense Frame Dimensions	900mm*500mm*2005mm
Built-in Dimensions	403mm*478mm*1600mm
Handle	ABS chrome plating
Slide Storage Capacity	129600 blocks/group



GROSSING STATION

- Ergonomically designed, fully functional, allowing you to stay away from pollution and poisoning.
- Offering safe, convenient & comfortable working environment.
- Materials: Quality stainless steel plate with 1.2mm thickness, water proof, rustproof, aci-& alkali-resistant.



SPECIFICATION PARAMETERS

Product Name	Grossing Station
Item No.	313020001
Model	PW900
Dimensions	1800mm*850mm*2200mm
Item No.	313020002
Model	PW901
Dimensions	1500mm*850mm*2200mm

CRYO-CABINET SF736

- Cryo-Cabinet is for biopsy specimen storage in pathology.
- Single-compressor freezing, double security, effectively prevents formaldehyde volatilization.
- Double-door with large capacity storage space, removable specimen storage trays, easy to operate.



SPECIFICATION PARAMETERS

Item No.	313000001
Product Name	Cryo-Cabinet
Model	SF736
Dimensions	1200mm*650mm*1910mm
Line-in(voltage)	AC220V-50Hz
Rating	240W
Temperature	0°C~12°C(Adjustable)
Capacity	1200L

EMBEDDING CASSETTES

EXCEL EMBEDDING CASSETTES



Item No.	Colour	
31201010017	White <input type="checkbox"/>	250 pcs/box, 8 Boxes/Carton Stripe Holes, with removable lid. P.O.M.Material; can be used on cassette printing machines.
31201020017	Blue <input type="checkbox"/>	
31201030017	Pink <input type="checkbox"/>	
31201040017	Purple <input type="checkbox"/>	
31201060017	Yellow <input type="checkbox"/>	
31201070017	Green <input type="checkbox"/>	



Item No.	Colour	
31201010008	White <input type="checkbox"/>	4*500 bottoms, 2*1000 lids/carton Stripe Holes, with detachable lid, Assembled. P.O.M. Material; can be used on cassette printing machines.
31201020008	Blue <input type="checkbox"/>	
31201030008	Pink <input type="checkbox"/>	
31201040008	Purple <input type="checkbox"/>	
31201060017	Yellow <input type="checkbox"/>	
31201070008	Green <input type="checkbox"/>	



Item No.	Colour	
31201010011	White <input type="checkbox"/>	250 pcs/box, 8 Boxes/Carton Square Holes, with removable lid. P.O.M.Material; can be used on cassette printing machines.
31201020011	Blue <input type="checkbox"/>	
31201030011	Pink <input type="checkbox"/>	
31201040011	Purple <input type="checkbox"/>	
31201060011	Yellow <input type="checkbox"/>	
31201070011	Green <input type="checkbox"/>	



Item No.	Colour	
31201010012	White <input type="checkbox"/>	250 pcs/box, 8 Boxes/Carton Fine checkered Holes, with detachable lid, Assembled. P.O.M. Material; can be used on cassette printing machines.
31201020012	Blue <input type="checkbox"/>	
31201030012	Pink <input type="checkbox"/>	
31201040012	Purple <input type="checkbox"/>	
31201060012	Yellow <input type="checkbox"/>	
31201070012	Green <input type="checkbox"/>	

Item No.	Colour	
31201010022	White <input type="checkbox"/>	4*500 bottoms, 2*1000 lids/carton Stripe Holes, with vertically rear lid. P.O.M.Material; can be used on cassette printing machines.
31201020020	Blue <input type="checkbox"/>	
31201030020	Pink <input type="checkbox"/>	
31201040020	Purple <input type="checkbox"/>	
31201060020	Yellow <input type="checkbox"/>	
31201070020	Green <input type="checkbox"/>	



Item No.	Colour	
31201010009	White <input type="checkbox"/>	4*500 bottoms, 2*1000 lids/carton Stripe Holes, with rear lid. P.O.M.Material; can be used on cassette printing machines.
31201020009	Blue <input type="checkbox"/>	
31201030009	Pink <input type="checkbox"/>	
31201040009	Purple <input type="checkbox"/>	
31201060009	Yellow <input type="checkbox"/>	
31201070009	Green <input type="checkbox"/>	



Item No.	Colour	
31201010007	White <input type="checkbox"/>	500 pcs/box, 8 Boxes/Carton Round Holes, without lid. P.O.M. Material; can be used on cassette printing machines.
31201020007	Blue <input type="checkbox"/>	
31201030007	Pink <input type="checkbox"/>	
31201040007	Purple <input type="checkbox"/>	
31201060007	Yellow <input type="checkbox"/>	
31201070007	Green <input type="checkbox"/>	



EMBEDDING CASSETTES

LASER EMBEDDING CASSETTES



Item No.	Colour	
31291010002	White <input type="checkbox"/>	4*500 bottoms, 2*1000 lids/carton Laser Embedding Cassettes, with rear lid. P.O.M.Material; can be used on Laser printing machines.
31291020002	Blue <input type="checkbox"/>	
31291030002	Pink <input type="checkbox"/>	
31291040002	Purple <input type="checkbox"/>	
31291060002	Yellow <input type="checkbox"/>	
31291070002	Green <input type="checkbox"/>	



Item No.	Colour	
31291010001	White <input type="checkbox"/>	4*500 bottoms, 2*1000 lids/carton Laser Embedding Cassettes, with detachable lid, Assembled. P.O.M.Material; can be used on Laser printing machines.
31291020001	Blue <input type="checkbox"/>	
31291030001	Pink <input type="checkbox"/>	
31291040001	Purple <input type="checkbox"/>	
31291060001	Yellow <input type="checkbox"/>	
31291070001	Green <input type="checkbox"/>	

PRO EMBEDDING CASSETTES



Item No.	Colour	
31201010014	White <input type="checkbox"/>	4*500 bottoms, 2*1000 lids/carton Stripe Holes, Unique design with sliding slot. can be used on cassette printing machines.
31201020014	Blue <input type="checkbox"/>	
31201030014	Pink <input type="checkbox"/>	
31201040014	Purple <input type="checkbox"/>	
31201060014	Yellow <input type="checkbox"/>	
31201070014	Green <input type="checkbox"/>	

Item No.	Colour	
31201010015	White <input type="checkbox"/>	250 pcs/box, 8 Boxes/Carton Gastroscopic tissue Embedding Cassettes, can be used on cassette printing machines.
31201020015	Blue <input type="checkbox"/>	
31201030015	Pink <input type="checkbox"/>	
31201040015	Purple <input type="checkbox"/>	
31201060015	Yellow <input type="checkbox"/>	
31201070015	Green <input type="checkbox"/>	



Item No.	Item No.	Colour	
31201010001	31201010002	White <input type="checkbox"/>	w/o lid: 4000 pcs/Carton lid: 2000 pcs/Carton Stripe Holes. Unique design with sliding slot, easy assembly & disassembly of the lids.
31201020001	31201020002	Blue <input type="checkbox"/>	
31201030001	31201030002	Pink <input type="checkbox"/>	
31201040001	31201040002	Purple <input type="checkbox"/>	
31201060001	31201060002	Yellow <input type="checkbox"/>	
31201070001	31201070002	Green <input type="checkbox"/>	



TRANSPARENT TUBULAR

Item No.	Item No.	Colour	
31201010003	31201010004	White <input type="checkbox"/>	w/o lid: 4500 pcs/Carton lid: 3000 pcs/Carton Stripe Holes. Working with cassette labeler; Quality materials of the tubes, deformation-resistant & smooth transport.
31201020003	31201020004	Blue <input type="checkbox"/>	
31201030003	31201030004	Pink <input type="checkbox"/>	
31201040003	31201040004	Purple <input type="checkbox"/>	
31201060003	31201060004	Yellow <input type="checkbox"/>	
31201070003	31201070004	Green <input type="checkbox"/>	



STAINLESS STEEL EMBEDDING MOLDS

- Quality stainless steel, durable;
- Four different specifications, depending design in sample trough, suitable for all types of specimens of different sizes.



STAINLESS LIDS



SPECIFICATION PARAMETERS

Item No.	31201000002
Dimensions	37mm*24mm*6mm
Packing	100pcs/box

Item No.	31201000003
Dimensions	30mm*24mm*6mm
Packing	100pcs/box

Item No.	31201000004
Dimensions	24mm*24mm*6mm
Packing	100pcs/box

Item No.	31201000005
Dimensions	10mm*10mm*5mm
Packing	100pcs/box

Item No.	31201000001
Packing	100pcs/box




DRYING PLATE

- Environmental protection high quality ABS material, strong and durable corrosion resistance.
- Identification label can be attached to the bottom.
- Ventilated and breathable, conducive to slide drying.
- A variety of colors to meet the needs of experiment.




- Dimensions: 308mm*194mm*11mm
- Pack of 20pcs, 60pcs/Carton

SPLIT-TYPE DRYING PLATE



Item No.	Colour
31205010001	White <input type="checkbox"/>
31205020001	Blue <input type="checkbox"/>
31205040001	Purple <input type="checkbox"/>
31205050001	Orange <input type="checkbox"/>
31205060001	Yellow <input type="checkbox"/>
31205070001	Green <input type="checkbox"/>
31205080001	Red <input type="checkbox"/>
31205000001	Quicko (mixed loading)

FLAT-TYPE DRYING PLATE



Item No.	Colour
31205010002	White <input type="checkbox"/>
31205020002	Blue <input type="checkbox"/>
31205040002	Purple <input type="checkbox"/>
31205050002	Orange <input type="checkbox"/>
31205060002	Yellow <input type="checkbox"/>
31205070002	Green <input type="checkbox"/>
31205080002	Red <input type="checkbox"/>
31205000002	Quicko (mixed loading)

PATHOLOGY GRADE PARAFFIN WAX

- Good toughness, offers ultra-thin and smooth section with no wrinkles.
- Easy dewaxing, ideal staining result.
- High purity, with no odor.
- Oil content: up to 0.5%.




SPECIFICATION PARAMETERS

Item No.	312020001
Melting Point	54~56°C
Colour	White
Packing	1kg/Pack; 8Packs/Carton

Item No.	312020002
Melting Point	57~60°C
Colour	White
Packing	1kg/Pack; 8Packs/Carton

DISSECTION BLADE HANDLE

- The overall stainless steel material for durability.
- For the sampling of tissue specimens in the pathology department.
- Blade fasten tightly, changing blade is safe and convenient.
- Suitable for a variety of blades to achieve blade reuse.



SPECIFICATION PARAMETERS

Item No.	Model
313020003	W80E
313020004	W80M
313020006	W80F
313020007	W80EM

MICROSCOPE GLASS SLIDES

STANDARD MICROSCOPE GLASS SLIDES



Item No.	312030001
Dimensions	75mm*25mm
Thickness	1.0mm~1.2mm
Packing	50Pcs/Box; 50Boxes/Carton

Flotation Glass Quality, Single-Head and Single-Side Frosting, 45° Angle; Polished Edge.



Item No.	312030002
Dimensions	75mm*25mm
Thickness	1.0mm~1.2mm
Packing	50Pcs/Box; 50Boxes/Carton

Flotation Glass Quality, Single-Head and Two-Sided Frosting, 45° Angle; Polished Edge.

EXCEL MICROSCOPE GLASS SLIDES



Item No.	312030003
Dimensions	75mm*25mm
Thickness	1.0mm~1.2mm
Packing	50Pcs/Box; 40Boxes/Carton

Super Grade Green Glass Material; Single Frosted; Ground Edges; 45° Corners.



Item No.	312030006
Dimensions	75mm*25mm
Thickness	1.0mm~1.2mm
Packing	50Pcs/Box; 40Boxes/Carton

Super Grade Green Glass Material; Double Frosted; Ground Edges; 45° Corners.

WEXIS-PRO MICROSCOPE GLASS SLIDES



Item No.	312030004
Dimensions	75mm*25mm
Thickness	1.0mm~1.2mm
Packing	50Pcs/Box; 40Boxes/Carton

Super White Glass Material; Single Painted; Ground Edges; 45° Corners.

ADHESION MICROSCOPE GLASS SLIDES



Item No.	312030005
Dimensions	75mm*25mm
Thickness	1.0mm~1.2mm
Packing	50Pcs/Box; 40Boxes/Carton

Super White Glass Class Material, Pre-cleaned, Specialized manufactured for IHC and Other Cytological and Histological tests. Positive Charged, Single Painted; Ground Edges; 45° Corners.

COVER GLASS

- Unique premium super white glass material, highly clear and transparent.



SPECIFICATION PARAMETERS

Item No.	Dimensions	Packing
312040001	24mm*24mm	100Pcs/Box; 200Boxes/Carton
312040002	24mm*32mm	100Pcs/Box; 200Boxes/Carton
312040003	24mm*40mm	100Pcs/Box; 200Boxes/Carton
312040003	24mm*50mm	100Pcs/Box; 200Boxes/Carton

FROZEN TISSUE EMBEDDING CENTER

- It is soluble in water and has no effect on dyeing.
- It supports the tissue effectively and reduces wrinkles and fragmentation when sectioned.



SPECIFICATION PARAMETERS

Item No.	Capacity
312000001	50ml/bottle
312000002	100ml/bottle
312000003	120ml/bottle
312000004	200ml/bottle

TISSUE FIXATIVE

- The reagent is safe, save time, labor friendly and reduce the air pollution of working environment effectively.
- Consistence and standardized, quality guaranteed.
- The reagent is per-packed in reagent bottles/containers for immediate use. No fur ther preparation required.
- Particularly applicable for immuno-histochemical staining of tissue. Better preservation of antigen of cells. Clearer and easier localization of antigen.
- Satisfies quality control requirement of immuno-histochemical technique.
- Avoid the cross-contamination and Malpractice.
- Reduce costs and improve efficacy.



SPECIFICATION PARAMETERS

Item No.	Capacity	Packing
311010001	6ml/bottle	500btl/Carton
311010004	15ml/bottle	500btl/Carton
311010005	20ml/bottle	500btl/Carton
311010006	30ml/bottle	320btl/Carton
311010007	50ml/bottle	200btl/Carton
311010008	100ml/bottle	180btl/Carton
311010009	250ml/bottle	50btl/Carton
311010010	500ml/bottle	20btl/Carton
311010011	1000ml/bottle	10btl/Carton
311010012	2500ml/bottle	8btl/Carton
311010013	5000ml/bottle	4btl/Carton
311010014	10L/pail	1pail/Carton
311010015	15L/pail	1pail/Carton
311010016	20L/pail	1pail/Carton

HEMATOXYLIN-EOSIN SOLUTION

- Hematoxylin solution is fuchsia, deposit-free, having a slightly pungent odour.
- Eosin solution is orange-red or purplish-red, deposit-free having a slightly pungent odour.
- A stained typical tissue observed under the microscope at room temperature, nuclei are stained blue, whereas the cytoplasm and extracellular matrix have varying degrees of pink staining. Well-fixed cells show considerable intramuscular detail.



SPECIFICATION PARAMETERS

Item No.	Description	Capacity	Packing
311020001	Lillie-Mayer's Hematoxylin Solution	500ml/bottle	6btl's/Carton
311020002	Lillie-Mayer's Hematoxylin Solution	2500ml/bottle	2btl's/Carton
311020003	Mayer's Hematoxylin Solution	500ml/bottle	6btl's/Carton
311020004	Mayer's Hematoxylin Solution	2500ml/bottle	2btl's/Carton
311020005	Harri's Hematoxylin Solution	500ml/bottle	6btl's/Carton
311020006	Harri's Hematoxylin Solution	2500ml/bottle	2btl's/Carton
311020007	Eosin Y, alcohol soluble	500ml/bottle	6btl's/Carton
311020008	Eosin Y, alcohol soluble	2500ml/bottle	2btl's/Carton
311020009	Eosin Y, water soluble	500ml/bottle	6btl's/Carton
311020010	Eosin Y, water soluble	2500ml/bottle	2btl's/Carton

SPECIAL STAIN SOLUTION



SPECIFICATION PARAMETERS

Item No.	Product Name	Specification	Intended Use
311020015	Collagen Fibers and Muscle Fibers Staining Kit(Van Gieson)	8ml/bottle*4	Van Gieson's Method for the detection of Collagen Fibers, muscle and cornified epithelium.
311020016	Collagen Fibers and Muscle Fibers Staining Kit(Masson)	8ml/bottle*5	Detection of Collagen Fibers in tissues such as skin, heart,etc
311020017	Reticular Fibers Staining Kit (Silver Hydroxide & Ammonia)	8ml/bottle*5	Detection of Reticular Fibers in tissues
311020018	Elastic Fibers Staining Kit(Orcein)	8ml/bottle*4	Detection of Elastic Fibers in tissues
311020019	Muscle Fibers Staining Kit(PTAH)	8ml/bottle*5	Detection of Muscle Fibers in tissues
311020020	Acid-fast Bacilli Staining Kit (Phenol Fuchsin)	8ml/bottle*3	Detection of Acid-fast bacilli in tissues
311020021	Fungi Staining Kit(Hexamine Silver)	8ml/bottle*12	Detection of Fungi in tissues
311020022	Amyloid Protein Staining Kit (Methanol Congo Red)	8ml/bottle*4	Detection of Amyloid Protein in tissues
311020023	Melanin Staining Kit(Ferrous Sulfate)	8ml/bottle*3+20bags	Detection of Melanin in tissues
311020024	Hemosiderin Staining Kit (Potassium Ferrocyanide)	8ml/bottle*3	Detection of Hemosiderin in tissues
311020025	Calcium Salt Staining Kit(Silver Nitrate)	8ml/bottle*4	Demonstrating deposits of calcium or calcium salt
311020026	Pituitary Cells Staining Kit(Periodic Acid-Colorless Magenta-Orange)	8ml/bottle*5	Detection of Pituitary Cells in tissues
311020027	Neutral Fat Staining Kit(Sudan III)	8ml/bottle*3	Detection of Neutral Fat in tissues
311020028	HBsAg Staining Kit(Aldehyde Fuchsin)	8ml/bottle*5	Detection of Hepatitis B surface antigen (HBsAg) in tissues
311020029	Helicobacter Pylori(HP) Staining Kit(BAMB)	8ml/bottle*1	Detection of Helicobacter Pylori (HP) in tissues
311020030	Lipofuscin Staining Kit(Aldehyde Fuchsin)	8ml/bottle*5	Detection of Lipofuscin in tissues
311020031	Mast Cell Staining Kit(Aldehyde Fuchsin)	8ml/bottle*3	Staining examination of histological mast cell
311020032	Myelin Sheath Staining Kit (Solochrome Cyanin)	8ml/bottle*4	Detection of Myelin Sheath in tissues
311020033	DNA Staining Kit(acid hydrolysis-Colorless Magenta)	8ml/bottle*4	Detection of Deoxyribonucleotides intissues
311020034	Glycogen Staining Kit(Periodic Acid- Colorless Magenta)	8ml/bottle*5	Staining glycosylated proteins in polyacrylamide gels usingthe periodic acid-Schiff(PAS)method
311020035	Mucosubstance Staining Kit (Alcian Blue(PH2.5))	8ml/bottle*2	Histological visualization of sulfated and carboxylated acid mucopolysaccharides and sulfated and carboxylated sialomucins(glycoproteins)
311020036	Mucosubstance Staining Kit (Alcian Blue(PH1.0))	8ml/bottle*2	Histological visualization of strongly sulfated mucosubstances
311020037	Mucosubstance Staining Kit (HID-AB)	8ml/bottle*3	Differentiation of sulfuric acid mucus and sialic acid mucus.
311020038	Basement Membrane Staining Kit (Hexamine Silver)	8ml/bottle*12	Demonstrating the basement membrane of the glomerulus in the kidney
311020039	Quick and Simple Acid Fast Staining Kit	8ml/bottle*6	Differentiate bacteria into acid fast group and non-acid fast groups
311020040	Carcinoid Staining Kit(Buered Silver Nitrate)	8ml/bottle*3+20bags	Detection of Carcinoid in tissues
311020041	Neutral Mucosubstance and Acid Mucosubstance Staining Kit(AB-PAS)	8ml/bottle*6	Detection of Neutral Mucosubstance and Acid Mucosubstance in tissues

HYALINE LIQUID&DEHYDRATED LIQUID



SPECIFICATION PARAMETERS

Item No.	Product Name	Description
311000013	Hyaline Liquid	/
311000014	Dehydrated Liquid	100% Ethanol
311000015	Dehydrated Liquid	95% Ethanol
311000018	Dehydrated Liquid	85% Ethanol
311000019	Dehydrated Liquid	75% Ethanol

Item No.	Capacity	Packing
311000013	500ml/bottle	20btl/Carton
311000014	500ml/bottle	20btl/Carton
311000015	500ml/bottle	20btl/Carton
311000018	500ml/bottle	20btl/Carton
311000019	500ml/bottle	20btl/Carton

OTHER REAGNT



SPECIFICATION PARAMETERS

Item No.	Product Name	Product Name	Capacity	Packing
311000001	Decolorization Solution	Decolorization Solution	500ml/bottle	6btl/Carton
311000003	Blue stain solution	Blue stain solution	500ml/bottle	6btl/Carton
311000004	Bouin's Solution	Bouin's Solution	500ml/bottle	6btl/Carton
311000005	Decalcifying Solution	Decalcifying Solution	500ml/bottle	6btl/Carton
311000006	Decalcifying Solution	Decalcifying Solution	2500ml/bottle	2btl/Carton
311000007	Decalcifying Solution	For Bone Marrow	500ml/bottle	6btl/Carton
311000008	Decalcifying Solution	For Bone Marrow	2500ml/bottle	2btl/Carton
311000009	Decalcifying Solution	8% Hydrochloric Acid	500ml/bottle	6btl/Carton
311000010	Decalcifying Solution	8% Hydrochloric Acid	2500ml/bottle	2btl/Carton
311000011	Decalcifying Solution	/	500ml/bottle	6btl/Carton
311020011	Papanicoiau Stain	Hematoxylin/og-6/EA-65	500ml/bottle	3btl/Carton
311020012	Papanicoiau Stain	Hematoxylin/og-6/EA-36	500ml/bottle	3btl/Carton
311020013	Papanicoiau Stain	Hematoxylin/og-6/EA-31	500ml/bottle	3btl/Carton
311020014	Papanicoiau Stain	Hematoxylin/og-6/EA-50	500ml/bottle	3btl/Carton
311010021	Frozen Section Fixative	For Frozen Section	50ml/bottle	200btl/Carton
311010022	Frozen Section Fixative	For Frozen Section	500ml/bottle	6btl/Carton